

# Certificate of performance



ISO 9001  
CERTIFIED

Article: 3 670 018 P

Vol. flasks, 25 ml, NS 14, amber glass, A, conf. cert.

Batch: 33-25


Serial:

Measuring Nr. 74164

Batch size:

141

Volume:	25,0000 ml
Tolerance:	± 0,0600 ml
Average value:	24,9921 ml
Deviation:	0,0073 ml
Precision (CV):	0,0291 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:

7

Date:

13.08.2025

Balance:

Sartorius -  
LE244S-OCE

DE-M