

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 670 025 USP

Vol. flask, 250 ml, NS 14, amber gl., USP, A, conf. cert.

Batch: 49-25

Serial:

Measuring Nr. 78081

Batch size:

278


Volume: 250,0000 ml

Tolerance: $\pm 0,1200$ ml

Average value: 249,9769 ml

Deviation: 0,0372 ml

Precision (CV): 0,0149 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

4

Date:

21.05.2026

Balance:

Sartorius - LA
1200 S

DE-M