

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 670 050

**Vol. flask, 500 ml, NS 19, amber glass, A, conf. cert.**


**Batch:** 32-21

**Serial:**

**Measuring Nr.** 48804

**Batch size:** 387

<b>Volume:</b>	500,000 ml
<b>Tolerance:</b>	$\pm 0,2500$ ml
<b>Average value:</b>	499,9926 ml
<b>Deviation:</b>	0,0284 ml
<b>Precision (CV):</b>	0,0057 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:** 9

**Date:** 13.09.2021

**Balance:** Sartorius -  
MSE14202S-  
OCE-D0

DE-M