

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 670 050

Vol. flask, 500 ml, NS 19, amber glass, A, conf. cert.


Batch: 16-26

Serial:

Measuring Nr. 77604

Batch size: 280

| | |
|------------------------|-----------------|
| Volume: | 500,000 ml |
| Tolerance: | $\pm 0,2500$ ml |
| Average value: | 500,0383 ml |
| Deviation: | 0,0416 ml |
| Precision (CV): | 0,0083 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator: 4

Date: 23.04.2026

Balance: Sartorius - LA
1200 S

DE-M