

Certificate of performance



Article: 3 670 101

Vol. flasks, 1000 ml, NS 29, amber glass, A, conf. cert.


Batch: 30-25

Serial:

Measuring Nr. 74099

Batch size: 90

Volume:	1000,000 ml
Tolerance:	± 0,6000 ml
Average value:	999,9421 ml
Deviation:	0,0866 ml
Precision (CV):	0,0087 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	12.08.2025
Balance:	Sartorius - MSE14202S- OCE-D0

