

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 510 011 BL

**Grad. pipette, 0,1:0,01 ml, blue grad., AS, conf. cert.**




**Batch:** 12-20

**Serial:**

**Measuring Nr.** 36740

**Batch size:** 142

<b>Volume:</b>	<b>0,1000 ml</b>	<b>0,0600 ml</b>	<b>0,0020 ml</b>
<b>Tolerance:</b>	$\pm 0,0060$ ml	$\pm 0,0060$ ml	$\pm 0,0060$ ml
<b>Average value:</b>	0,1009 ml	0,0601 ml	0,0188 ml
<b>Deviation:</b>	0,0020 ml	0,0013 ml	0,0007 ml
<b>Precision (CV):</b>	1,9572 %	2,2262 %	36,5504 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	4
<b>Date:</b>	24.03.2020
<b>Balance:</b>	Sartorius - MSE224S-0CE -DU

