

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 510 011

Grad. pipettes, 0,1:0,01 ml, AS, conf. cert.




Batch: 06-25

Serial:

Measuring Nr. 71899

Batch size: 280

Volume:	0,1000 ml	0,0600 ml	0,0200 ml
Tolerance:	$\pm 0,0060$ ml	$\pm 0,0060$ ml	$\pm 0,0060$ ml
Average value:	0,0977 ml	0,0587 ml	0,0182 ml
Deviation:	0,0013 ml	0,0008 ml	0,0004 ml
Precision (CV):	1,3393 %	1,3355 %	1,7858 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

4

Date:

12.02.2025

Balance:

Sartorius -
MSE224S-0CE
-DU

DE-M