

# Certificate of performance



**Article:** 3 510 051 BL

**Grad. pipettes, 0,5:0,01 ml, blue grad., AS, conf. cert.**




**Batch:** 24-23

**Serial:**

**Measuring Nr.** 63439

**Batch size:** 280

Volume:	0,5000 ml	0,3000 ml	0,1000 ml
<b>Tolerance:</b>	± 0,0060 ml	± 0,0060 ml	± 0,0060 ml
<b>Average value:</b>	0,4991 ml	0,3007 ml	0,0997 ml
<b>Deviation:</b>	0,0019 ml	0,0016 ml	0,0008 ml
<b>Precision (CV):</b>	0,3811 %	0,5309 %	0,7965 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	7
<b>Date:</b>	21.06.2023
<b>Balance:</b>	Sartorius - LE244S-OCE

