

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 510 101 BL

**Grad. pipettes, 1:0,01 ml, blue grad., AS, conf. cert.**




**Batch:** 25-22

**Serial:**

**Measuring Nr.** 56625

**Batch size:** 500

Volume:	1,0000 ml	0,6000 ml	0,2000 ml
Tolerance:	± 0,0070 ml	± 0,0070 ml	± 0,0070 ml
Average value:	1,0001 ml	0,6011 ml	0,2018 ml
Deviation:	0,0016 ml	0,0016 ml	0,0011 ml
Precision (CV):	0,1646 %	0,2601 %	0,5730 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:** 7

**Date:** 21.06.2022

**Balance:** Sartorius -  
LE244S-OCE

DE-M