

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 510 101

Grad. pipettes, 1:0,01 ml, AS, conf. cert.

Batch: 33-25




Serial:

Measuring Nr. 75043

Batch size:

3200

Volume:	1,0000 ml	0,6000 ml	0,2000 ml
Tolerance:	± 0,0070 ml	± 0,0070 ml	± 0,0070 ml
Average value:	0,9971 ml	0,6018 ml	0,2011 ml
Deviation:	0,0017 ml	0,0022 ml	0,0010 ml
Precision (CV):	0,1718 %	0,3635 %	0,4798 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

4

Date:

27.10.2025

Balance:

Sartorius -
MSE224S-0CE
-DU

