

Certificate of performance



Article: 3 510 202

Grad. pipettes, 2:0,02 ml, AS, conf. cert.




Batch: 18-18

Serial:

Measuring Nr. 20506

Batch size: 8000

Volume:	2,0000 ml	1,2000 ml	0,4000 ml
Tolerance:	± 0,0100 ml	± 0,0100 ml	± 0,0100 ml
Average value:	1,9998 ml	1,1979 ml	0,4003 ml
Deviation:	0,0027 ml	0,0022 ml	0,0015 ml
Precision (CV):	0,1363 %	0,1822 %	0,3651 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	9
Date:	16.05.2018
Balance:	Sartorius - LE244S-OCE

