

Certificate of performance



Article: 3 510 210 BL

Grad. pipettes, 2:0,1 ml, blue grad., AS, conf. cert.




Batch: 10-24

Serial:

Measuring Nr. 66336

Batch size: 280

| Volume: | 2,0000 ml | 1,2000 ml | 0,4000 ml |
|------------------------|-------------|-------------|-------------|
| Tolerance: | ± 0,0100 ml | ± 0,0100 ml | ± 0,0100 ml |
| Average value: | 1,9969 ml | 1,2007 ml | 0,4010 ml |
| Deviation: | 0,0023 ml | 0,0032 ml | 0,0016 ml |
| Precision (CV): | 0,1140 % | 0,2641 % | 0,4065 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|------------------|---------------------------|
| Operator: | 7 |
| Date: | 12.03.2024 |
| Balance: | Sartorius - LE244S-OCE |

