

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 510 505

Graduated pipettes, 5:0,05 ml, AS, conf. cert.




Batch: 05-18

Serial:

Measuring Nr. 18279

Batch size: 3200

Volume:	5,0000 ml	3,0000 ml	1,0000 ml
Tolerance:	$\pm 0,0300$ ml	$\pm 0,0300$ ml	$\pm 0,0300$ ml
Average value:	4,9939 ml	2,9916 ml	0,9995 ml
Deviation:	0,0066 ml	0,0041 ml	0,0028 ml
Precision (CV):	0,1324 %	0,1361 %	0,2832 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

9

Date:

21.02.2018

Balance:

Sartorius -
LE244S-OCE

DE-M