

Certificate of performance



Article: 3 510 505

Graduated pipettes, 5:0,05 ml, AS, conf. cert.




Batch: 47-18

Serial:

Measuring Nr. 25877

Batch size: 2400

Volume:	5,0000 ml	3,0000 ml	1,0000 ml
Tolerance:	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
Average value:	4,9908 ml	2,9979 ml	1,0008 ml
Deviation:	0,0056 ml	0,0059 ml	0,0034 ml
Precision (CV):	0,1113 %	0,1955 %	0,3382 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	9
Date:	26.11.2018
Balance:	Sartorius - LE244S-OCE

