

Certificate of performance



Article: 3 510 505

Graduated pipettes, 5:0,05 ml, AS, conf. cert.




Batch: 12-19

Serial:

Measuring Nr. 28911

Batch size: 3200

Volume:	5,0000 ml	3,0000 ml	1,0000 ml
Tolerance:	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
Average value:	4,9889 ml	2,9983 ml	1,0003 ml
Deviation:	0,0068 ml	0,0065 ml	0,0050 ml
Precision (CV):	0,1368 %	0,2164 %	0,4999 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	9
Date:	02.04.2019
Balance:	Sartorius - LE244S-OCE

