

Certificate of performance



Article: 3 510 710 BL

Grad. pipettes, 10:0,1 ml, blue grad., AS, conf. cert.




Batch: 24-24

Serial:

Measuring Nr. 67770

Batch size: 1197

Volume:	10,0000 ml	6,0000 ml	2,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	9,9987 ml	5,9929 ml	2,0012 ml
Deviation:	0,0076 ml	0,0099 ml	0,0052 ml
Precision (CV):	0,0760 %	0,1657 %	0,2623 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	24.06.2024
Balance:	Sartorius - LE244S-OCE

