

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 510 710 USP

Graduated pipettes, 10:0,1 ml, USP, AS, conf.cert.




Batch: 25-23

Serial:

Measuring Nr. 63600

Batch size: 280

Volume:	10,0000 ml	6,0000 ml	2,0000 ml
Tolerance:	$\pm 0,0300$ ml	$\pm 0,0300$ ml	$\pm 0,0300$ ml
Average value:	9,9990 ml	5,9892 ml	1,9992 ml
Deviation:	0,0082 ml	0,0046 ml	0,0033 ml
Precision (CV):	0,0815 %	0,0762 %	0,1660 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

31.10.2024

Balance:

Sartorius -
LE244S-OCE

DE-M