

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 510 710

Graduated pipettes, 10:0,1 ml, AS, conf. cert.




Batch: 02-24

Serial:

Measuring Nr. 65865

Batch size: 4765

Volume:	10,000 ml	6,000 ml	2,000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	9,9948 ml	6,0015 ml	2,0022 ml
Deviation:	0,0090 ml	0,0078 ml	0,0048 ml
Precision (CV):	0,0900 %	0,1308 %	0,2407 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

01.02.2024

Balance:

Sartorius -
LE244S-OCE

DE-M