

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 510 710

**Graduated pipettes, 10:0,1 ml, AS, conf. cert.**




**Batch:** 43-25

**Serial:**

**Measuring Nr.** 76133

**Batch size:** 4800

Volume:	10,0000 ml	6,0000 ml	2,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	9,9918 ml	5,9973 ml	2,0003 ml
Deviation:	0,0168 ml	0,0165 ml	0,0058 ml
Precision (CV):	0,1682 %	0,2743 %	0,2918 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	4
Date:	13.01.2026
Balance:	Sartorius - MSE224S-0CE -DU

