

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 510 810 BL

Grad. pipettes, 20:01, blue grad., AS, conf. cert




Batch: 25-26

Serial:

Measuring Nr. 78672

Batch size: 280

Volume:	20,0000 ml	12,0000 ml	4,0000 ml
Tolerance:	$\pm 0,1000 \text{ ml}$	$\pm 0,1000 \text{ ml}$	$\pm 0,1000 \text{ ml}$
Average value:	19,9739 ml	11,9956 ml	4,0001 ml
Deviation:	0,0153 ml	0,0160 ml	0,0069 ml
Precision (CV):	0,0765 %	0,1334 %	0,1726 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

4

Date:

24.06.2026

Balance:

Sartorius -
MSE224S-0CE
-DU

DE-M