

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 510 910 BL

Grad. pipettes, 25:0,1 ml, blue grad., AS, conf. cert.




Batch: 12-26

Serial:

Measuring Nr. 77973

Batch size: 280

Volume:	25,0000 ml	15,0000 ml	5,0000 ml
Tolerance:	$\pm 0,1000$ ml	$\pm 0,1000$ ml	$\pm 0,1000$ ml
Average value:	25,0462 ml	15,0449 ml	5,0221 ml
Deviation:	0,0210 ml	0,0184 ml	0,0105 ml
Precision (CV):	0,0841 %	0,1226 %	0,2090 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

06.05.2026

Balance:

Sartorius -
LE244S-OCE

DE-M