

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 510 910

**Grad. pipettes, 25:0,1 ml, AS, conf. cert.**




**Batch:** 38-24

**Serial:**

**Measuring Nr.** 70019

**Batch size:** 500

Volume:	25,0000 ml	15,0000 ml	5,0000 ml
<b>Tolerance:</b>	± 0,1000 ml	± 0,1000 ml	± 0,1000 ml
<b>Average value:</b>	24,9769 ml	15,0054 ml	5,0036 ml
<b>Deviation:</b>	0,0212 ml	0,0148 ml	0,0089 ml
<b>Precision (CV):</b>	0,0847 %	0,0985 %	0,1779 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	4
<b>Date:</b>	08.11.2024
<b>Balance:</b>	Sartorius - MSE224S-0CE -DU

