

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 510 910

Grad. pipettes, 25:0,1 ml, AS, conf. cert.




Batch: 12-26

Serial:

Measuring Nr. 78007

Batch size: 499

Volume:	25,0000 ml	15,0000 ml	5,0000 ml
Tolerance:	$\pm 0,1000$ ml	$\pm 0,1000$ ml	$\pm 0,1000$ ml
Average value:	25,0325 ml	15,0291 ml	5,0167 ml
Deviation:	0,0139 ml	0,0155 ml	0,0096 ml
Precision (CV):	0,0557 %	0,1031 %	0,1910 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

07.05.2026

Balance:

Sartorius -
LE244S-OCE

DE-M