

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 515 051

**Graduated pipettes, 0,5:0,01 ml, AS, conf.cert.**




**Batch:** 45-19

**Serial:**

**Measuring Nr.** 34404

**Batch size:** 150

Volume:	0,5000 ml	0,3000 ml	0,1000 ml
Tolerance:	± 0,0060 ml	± 0,0060 ml	± 0,0060 ml
Average value:	0,5005 ml	0,3011 ml	0,1002 ml
Deviation:	0,0017 ml	0,0005 ml	0,0003 ml
Precision (CV):	0,3340 %	0,1602 %	0,3412 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	7
Date:	12.11.2019
Balance:	Sartorius - LE244S-OCE

DE-M