

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 515 051

Graduated pipettes, 0,5:0,01 ml, AS, conf.cert.




Batch: 17-26

Serial:

Measuring Nr. 77952

Batch size: 150

Volume:	0,5000 ml	0,3000 ml	0,1000 ml
Tolerance:	$\pm 0,0060$ ml	$\pm 0,0060$ ml	$\pm 0,0060$ ml
Average value:	0,4970 ml	0,2961 ml	0,0988 ml
Deviation:	0,0010 ml	0,0007 ml	0,0006 ml
Precision (CV):	0,2093 %	0,2450 %	0,5740 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	05.05.2026
Balance:	Sartorius - MSE224S-0CE -DU

DE-M