

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 515 110

Grad. pipette, 1:0,1 ml, AS, conf. cert.




Batch: 41-25

Serial:

Measuring Nr. 75200

Batch size: 150

| Volume: | 1,0000 ml | 0,6000 ml | 0,2000 ml |
|------------------------|-------------|-------------|-------------|
| Tolerance: | ± 0,0070 ml | ± 0,0070 ml | ± 0,0070 ml |
| Average value: | 0,9983 ml | 0,5999 ml | 0,1989 ml |
| Deviation: | 0,0020 ml | 0,0021 ml | 0,0009 ml |
| Precision (CV): | 0,2020 % | 0,3472 % | 0,4468 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|------------------|---------------------------|
| Operator: | 7 |
| Date: | 05.11.2025 |
| Balance: | Sartorius - LE244S-OCE |

