

# Certificate of performance



**Article:** 3 515 505

**Grad. pipette, 5:0,05 ml, AS, conf. cert.**




**Batch:** 18-19

**Serial:**

**Measuring Nr.** 29873

**Batch size:** 500

Volume:	5,0000 ml	3,0000 ml	1,0000 ml
Tolerance:	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
Average value:	5,0022 ml	3,0026 ml	1,0036 ml
Deviation:	0,0065 ml	0,0053 ml	0,0041 ml
Precision (CV):	0,1294 %	0,1765 %	0,4086 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	7
Date:	13.05.2019
Balance:	Sartorius - LE244S-OCE

