

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 515 505

**Grad. pipette, 5:0,05 ml, AS, conf. cert.**




**Batch:** 27-23

**Serial:**

**Measuring Nr.** 63820

**Batch size:** 500

Volume:	5,0000 ml	3,0000 ml	1,0000 ml
<b>Tolerance:</b>	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
<b>Average value:</b>	4,9901 ml	2,9947 ml	1,0001 ml
<b>Deviation:</b>	0,0073 ml	0,0058 ml	0,0038 ml
<b>Precision (CV):</b>	0,1460 %	0,1929 %	0,3811 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	7
<b>Date:</b>	19.07.2023
<b>Balance:</b>	Sartorius - LE244S-OCE

