

# Certificate of performance



**Article:** 3 515 710

**Grad. pipettes, 10:0,1 ml, AS, conf. cert.**




**Batch:** 26-20

**Serial:**

**Measuring Nr.** 38815

**Batch size:** 500

Volume:	10,000 ml	6,0000 ml	2,0000 ml
<b>Tolerance:</b>	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
<b>Average value:</b>	9,9847 ml	5,9931 ml	2,0012 ml
<b>Deviation:</b>	0,0112 ml	0,0057 ml	0,0045 ml
<b>Precision (CV):</b>	0,1115 %	0,0951 %	0,2250 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	7
<b>Date:</b>	23.06.2020
<b>Balance:</b>	Sartorius - MSE224S-0CE -DU

