

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 515 710

Grad. pipettes, 10:0,1 ml, AS, conf. cert.




Batch: 05-24

Serial:

Measuring Nr. 65507

Batch size: 500

Volume:	10,000 ml	6,0000 ml	2,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	9,9973 ml	5,9976 ml	2,0014 ml
Deviation:	0,0101 ml	0,0092 ml	0,0046 ml
Precision (CV):	0,1011 %	0,1536 %	0,2323 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

07.02.2024

Balance:

Sartorius -
LE244S-OCE

DE-M