

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 515 710

**Grad. pipettes, 10:0,1 ml, AS, conf. cert.**

**Batch:** 20-26




**Serial:**

**Measuring Nr.** 78424

**Batch size:**

1200

|                        |                   |                  |                  |
|------------------------|-------------------|------------------|------------------|
| <b>Volume:</b>         | <b>10,0000 ml</b> | <b>6,0000 ml</b> | <b>2,0000 ml</b> |
| <b>Tolerance:</b>      | $\pm 0,0500$ ml   | $\pm 0,0500$ ml  | $\pm 0,0500$ ml  |
| <b>Average value:</b>  | 9,9886 ml         | 5,9932 ml        | 2,0014 ml        |
| <b>Deviation:</b>      | 0,0151 ml         | 0,0087 ml        | 0,0040 ml        |
| <b>Precision (CV):</b> | 0,1512 %          | 0,1446 %         | 0,1993 %         |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

7

**Date:**

15.06.2026

**Balance:**

Sartorius -  
LE244S-OCE

DE-M