

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 515 810

Grad. pipettes, 20:0,1 ml, AS, conf. cert.




Batch: 33-24

Serial:

Measuring Nr. 68850

Batch size: 150

Volume:	20,000 ml	12,000 ml	4,000 ml
Tolerance:	± 0,1000 ml	± 0,1000 ml	± 0,1000 ml
Average value:	19,9788 ml	11,9843 ml	3,9942 ml
Deviation:	0,0173 ml	0,0133 ml	0,0073 ml
Precision (CV):	0,0865 %	0,1111 %	0,1815 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	23.08.2024
Balance:	Sartorius - MSE224S-0CE -DU

