

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 515 910

**Graduated pipettes class AS 25 ml**




**Batch:** 33-25

**Serial:**

**Measuring Nr.** 74102

**Batch size:** 279

Volume:	25,000 ml	15,000 ml	5,000 ml
<b>Tolerance:</b>	± 0,1000 ml	± 0,1000 ml	± 0,1000 ml
<b>Average value:</b>	25,0470 ml	15,0305 ml	5,0239 ml
<b>Deviation:</b>	0,0207 ml	0,0110 ml	0,0130 ml
<b>Precision (CV):</b>	0,0828 %	0,0734 %	0,2605 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	7
<b>Date:</b>	13.08.2025
<b>Balance:</b>	Sartorius - LE244S-OCE

