

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 517 011

Pipette f. encym. analysis., 0,1:0,001 ml, AS, conf. cert.




Batch: 32-23

Serial:

Measuring Nr.: 64034

Batch size: 500

Volume:	0,1000 ml	0,0600 ml	0,0200 ml
Tolerance:	± 0,0030 ml	± 0,0030 ml	± 0,0030 ml
Average value:	0,0994 ml	0,0589 ml	0,0187 ml
Deviation:	0,0008 ml	0,0005 ml	0,0005 ml
Precision (CV):	0,8003 %	0,9125 %	2,6285 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to DIN 12699

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	15.08.2023
Balance:	Sartorius - LE244S-OCE

