

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 517 011

**Pipette f. encym. analysis., 0,1:0,001 ml, AS, conf. cert.**




**Batch:** 06-25

**Serial:**

**Measuring Nr.:** 71878

**Batch size:** 280

Volume:	0,1000 ml	0,0600 ml	0,0200 ml
Tolerance:	± 0,0030 ml	± 0,0030 ml	± 0,0030 ml
Average value:	0,1004 ml	0,0595 ml	0,0186 ml
Deviation:	0,0009 ml	0,0006 ml	0,0004 ml
Precision (CV):	0,9320 %	1,0737 %	2,1313 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to DIN 12699

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	4
Date:	03.02.2025
Balance:	Sartorius - MSE224S-0CE -DU

