

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 517 021

Pipette f. encym. analysis., 0,2:0,001 ml, AS, conf. cert.




Batch: 17-20

Serial: -

Measuring Nr. 38027

Batch size: 90

Volume:	0,2000 ml	0,1200 ml	0,0400 ml
Tolerance:	± 0,0030 ml	± 0,0030 ml	± 0,0030 ml
Average value:	0,1997 ml	0,1185 ml	0,0388 ml
Deviation:	0,0007 ml	0,0005 ml	0,0005 ml
Precision (CV):	0,3263 %	0,4001 %	1,2792 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to DIN 12699

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	18.05.2020
Balance:	Sartorius - MSE224S-0CE -DU

