

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 517 101 BL

Pipettes f. encym. analysis, 1:0,01 ml, AS, conf. cert.




Batch: 36-20

Serial:

Measuring Nr. 39986

Batch size: 90

Volume:	1,0000 ml	0,6000 ml	0,2000 ml
Tolerance:	± 0,0060 ml	± 0,0060 ml	± 0,0060 ml
Average value:	1,0002 ml	0,6003 ml	0,2006 ml
Deviation:	0,0023 ml	0,0015 ml	0,0012 ml
Precision (CV):	0,2344 %	0,2484 %	0,6098 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to DIN 12699

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	02.09.2020
Balance:	Sartorius - LE244S-OCE

DE-M