

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 517 101

**Pipettes f. encym. analysis, 1:0,01 ml, AS, conf. cert.**

**Batch:** 22-20




**Serial:**

**Measuring Nr.** 38345

**Batch size:**

249

Volume:	1,0000 ml	0,6000 ml	0,2000 ml
Tolerance:	± 0,0060 ml	± 0,0060 ml	± 0,0060 ml
Average value:	1,0011 ml	0,6008 ml	0,2004 ml
Deviation:	0,0022 ml	0,0015 ml	0,0010 ml
Precision (CV):	0,2179 %	0,2501 %	0,4890 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to DIN 12699

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	7
Date:	29.05.2020
Balance:	Sartorius - MSE224S-0CE -DU

