

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 517 101

**Pipettes f. encym. analysis, 1:0,01 ml, AS, conf. cert.**




**Batch:** 26-26

**Serial:**

**Measuring Nr.** 78722

**Batch size:** 280

|                        |                  |                  |                  |
|------------------------|------------------|------------------|------------------|
| <b>Volume:</b>         | <b>1,0000 ml</b> | <b>0,6000 ml</b> | <b>0,2000 ml</b> |
| <b>Tolerance:</b>      | $\pm 0,0060$ ml  | $\pm 0,0060$ ml  | $\pm 0,0060$ ml  |
| <b>Average value:</b>  | 1,0015 ml        | 0,6026 ml        | 0,2012 ml        |
| <b>Deviation:</b>      | 0,0018 ml        | 0,0013 ml        | 0,0008 ml        |
| <b>Precision (CV):</b> | 0,1834 %         | 0,2174 %         | 0,4051 %         |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to DIN 12699

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

7

**Date:**

29.06.2026

**Balance:**

Sartorius -  
LE244S-OCE

DE-M