

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 517 505 BL

Pipettes f. encym. analysis, 5:0,05 ml, AS, conf. cert.




Batch: 26-18

Serial:

Measuring Nr.: 21863

Batch size: 90

Volume:	5,0000 ml	3,0000 ml	1,0000 ml
Tolerance:	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
Average value:	5,0009 ml	2,9988 ml	1,0009 ml
Deviation:	0,0048 ml	0,0031 ml	0,0016 ml
Precision (CV):	0,0963 %	0,1024 %	0,1602 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to DIN 12699

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator: 44

Date: 02.07.2018

Balance: Sartorius -
LE244S-OCE

DE-M