

# Certificate of performance



**Article:** 3 517 505

**Pipettes for enzym. analysis, 5:0,05 ml, AS, conf. cert.**




**Batch:** 50-23

**Serial:**

**Measuring Nr.** 65493

**Batch size:** 100

Volume:	5,0000 ml	3,0000 ml	1,0000 ml
<b>Tolerance:</b>	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
<b>Average value:</b>	5,0016 ml	2,9976 ml	1,0005 ml
<b>Deviation:</b>	0,0066 ml	0,0049 ml	0,0019 ml
<b>Precision (CV):</b>	0,1310 %	0,1644 %	0,1889 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to DIN 12699

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	7
<b>Date:</b>	18.12.2023
<b>Balance:</b>	Sartorius - LE244S-OCE

