

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 535 010 USP

Volumetric pipettes, 1 ml, USP, AS, conf. cert.


Batch: 35-22

Serial:

Measuring Nr. 59371

Batch size: 279

| | |
|------------------------|-------------|
| Volume: | 1,0000 ml |
| Tolerance: | ± 0,0060 ml |
| Average value: | 0,9994 ml |
| Deviation: | 0,0013 ml |
| Precision (CV): | 0,1279 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 648

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|------------------|-----------------------------------|
| Operator: | 4 |
| Date: | 26.09.2022 |
| Balance: | Sartorius - MSE224S-0CE -DU |

DE-M