

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 535 070 USP

Volumetric pipette, 7 ml, USP, AS, conf. cert

Batch: 24-25


Serial:

Measuring Nr. 74034

Batch size:

146

| | |
|------------------------|-------------|
| Volume: | 7,0000 ml |
| Tolerance: | ± 0,0100 ml |
| Average value: | 6,9991 ml |
| Deviation: | 0,0030 ml |
| Precision (CV): | 0,0436 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 648

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|------------------|---------------------------|
| Operator: | 7 |
| Date: | 31.07.2025 |
| Balance: | Sartorius - LE244S-OCE |

DE-M