

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 535 400 BL

Volumetric pipettes, 40 ml, blue grad., AS, conf. cert.


Batch: 04-08

Serial:

Measuring Nr. 17755

Batch size: 180

| | |
|------------------------|------------------|
| Volume: | 40,000 ml |
| Tolerance: | $\pm 0,0500$ ml |
| Average value: | 39,9991 ml |
| Deviation: | 0,0088 ml |
| Precision (CV): | 0,0220 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 648

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator: 9

Date: 31.01.2018

Balance: Sartorius -
LE244S-OCE

DE-M