

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 555 005

Volumetric pipettes, 0,5 ml, 2 marks, AS, conf. cert.


Batch: 48-23

Serial:

Measuring Nr. 65425

Batch size: 474

| | |
|------------------------|------------------|
| Volume: | 0,5000 ml |
| Tolerance: | $\pm 0,0050$ ml |
| Average value: | 0,5011 ml |
| Deviation: | 0,0009 ml |
| Precision (CV): | 0,1729 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 648

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator: 7

Date: 05.12.2023

Balance: Sartorius -
LE244S-OCE

DE-M