

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 555 020

Volumetric pipettes, 2 ml, with 2 marks, AS, conf. cert.

Batch: 04-18


Serial:

Measuring Nr. 17474

Batch size:

1

| | |
|------------------------|-----------------|
| Volume: | 2,0000 ml |
| Tolerance: | $\pm 0,0100$ ml |
| Measured value: | 2,0001 ml |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 648

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|------------------|---------------------------|
| Operator: | 9 |
| Date: | 24.01.2018 |
| Balance: | Sartorius - LE244S-OCE |

DE-M