

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 555 300

Volumetric pipettes, 30 ml, 2 marks, AS, conf.cert.

Batch: 03-25

Serial:

Measuring Nr. 71329

Batch size:

150


Volume: 30,0000 ml

Tolerance: $\pm 0,0300$ ml

Average value: 29,9982 ml

Deviation: 0,0059 ml

Precision (CV): 0,0197 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 648

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

4

Date:

14.01.2025

Balance:

Sartorius -
MSE224S-0CE
-DU

DE-M