

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 555 500

Vol. pipettes, 50 ml, with 2 marks, AS, conf.cert.

Batch: 46-17


Serial:

Measuring Nr. 16530

Batch size:

279

Volume:	50,0000 ml
Tolerance:	$\pm 0,0500$ ml
Average value:	49,9983 ml
Deviation:	0,0084 ml
Precision (CV):	0,0169 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 648

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

9

Date:

18.12.2017

Balance:

Sartorius -
LE244S-OCE

DE-M