

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 555 500

Vol. pipettes, 50 ml, with 2 marks, AS, conf.cert.

Batch: 03-20


Serial:

Measuring Nr. 35844

Batch size:

199

Volume:	50,0000 ml
Tolerance:	$\pm 0,0500$ ml
Average value:	49,9943 ml
Deviation:	0,0099 ml
Precision (CV):	0,0197 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 648

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

04.02.2020

Balance:

Sartorius -
MSE224S-0CE
-DU

DE-M